


Ammonia A023

Example 1 - Range: 0.001 - 0.05 mg/L




Ammonia A023-A01

Range: 0.001 - 0.05 mg/L (as N)

MDL=0.0004 mg/L (as N)

Analysis rate: 45 samples/hour

Example 2 - Range: 0.05 - 25 mg/L




Ammonia A023-AMT

Range: 0.05 - 25 mg/L (as N)

MDL=0.007 mg/L (as N)

Analysis rate: 48 samples/hour

Example 3 - Range: 0.2 - 80 mg/L



Ammonia A023-BMT

Range: 0.2 - 80 mg/L (as N)

MDL=0.05 mg/L (as N)

Analysis rate: 60 samples/hour