

ESSENTIALS FOR THE WINE LABORATORY

Astoria-Pacific is positioning itself as a one-stop analytical partner for wineries and beverage producers, expanding our portfolio to include essential laboratory and production measurement tools. By combining practical instrumentation with technical expertise and responsive service, we help producers protect quality, optimize processes, and ensure confidence from vineyard to bottle.



The rAPID-T System automates common tests: Acetic Acid, Glucose+ Fructose, L-Malic Acid, NOPA, Ammonia, Total Sulfites, etc.
Free Sulfite available with Optional Sample Prep system.

The Astoria2 Analyzer automates the simultaneous analysis of Free and Total Sulfites at 120 tests/hr. Includes a sampler with capacity of 180 samples.



Astoria-Pacific has partnered with Alpine Scientific to service and support Check Stab Tartaric Stability, Protein, and Filterability systems from Delta Acque. Contact Alpine and/or us for availability in your area.



Titra EVO ACID Automatic titrator pH, TA, SO₂, & VA after extraction system from Dujardin-Salleron. *Optional 15-sample autosampler available for increased throughput and efficiency.*



DS Titra. Three different models from Dujardin-Salleron: (1) Sulfi, (2) Iodo, (3) Titra. Select the best instrument for the analysis of Free and Total Sulfites, or pH, TA and Acetic Acid (after extraction).

The electric ebulliometer from Dujardin-Salleron accurately determines the alcohol content (0 – 17%) of dry wines in ~ 7 min. ***Dujardin-Salleron's traditional ebulliometer is also available.*



Free Sulfite Sample Prep Station. Run FSO₂ analysis with confidence on your rAPID-T, spec, and/or your current automated analyzer. Long-lasting reagents. Uses the same method as the best reference labs like myEnologist and Baker Wine Lab.



Distillation, Extraction, Dissolved CO₂ systems from Dujardin-Salleron are available. Dujardin-Salleron conceives, manufactures and distributes analysis instruments and measurement devices for oenology and spirits.



Astoria-Pacific offers a variety of wine analysis reagent kits and consumables for its rAPID-T and Astoria Analyzers. We supply reliable and affordable solutions for your laboratory.

Astoria-Pacific is excited to offer products from Dujardin-Salleron. They have been providing wine laboratories reference equipment for over 160 years.



ASTORIA • PACIFIC